BIO-MAPPER

Product Data Report

Product Name: C-CRP Antigen

Catalog: APR130201

Host Animal: Ecoli, Recombiant Antigen

Fragment: Canine CRP gene fragment

Lot: 210701

Mfg Date: 2021-07-01

Exp Date: 2024-06-30

Con: 1.5 mg/ml

Purity: >95%, by SDS-PAGE, Affinity chromatography

Buffer: 50mM Tris-Cl pH8.0

Preservative: 0.1% Proclin300

Store: Store at -20°C±5°C.Repeated opening and closing of the bottle should be avoided. Antigens are stable until the

expiry date in an unopened vial.

Application: Later Flow/Enzyme Immunoassay/CMIA/CLIA

Usage: For application positive control. Each laboratory should determine an optimum working concentration for use in particular application.

Validity: 3 years.

Remark: This product is intended for research and manufacturing uses only. It is a diagnostic device. The user assumes all responsibility for care, custod and control of the material, including its disposal, in accordance with all regulations.

