## **BIO-MAPPER**

## **Product Data Report**

Product Name: HIV gp36 Antigen

Catalog: BMEHIV121

Host Animal: Ecoli ,Recombinant Antigen

Fragment: HIV M type gp41 gene fragment

**Lot:** 210701

Mfg Date: 2021-07-01

**Exp Date:** 2024-06-30

**Con:** 1.5 mg/ml

**Purity:** >95%,by SDS-PAGE,Affinity chromatography

**Buffer:** 50mM Tris-Cl pH10.0

**Preservative:** 0.1% Proclin300

Store: Store at 2-8 °C.Repeated opening and closing of the bottle should be avoided. Antigens are stable until the expi-

ry date in an unopened vial.

**Application:** Enzyme Immunoassay

Matched Pair: BMEHIV101/BMEHIV114

**Usage:** For application in indirect method for conjugation, pair with BMEHIV101 (HIV I+II fusion antigen for capture) and BMEHIV114 (HIV I-HRP conjugation) for HIV Antibody tesing. Each laboratory should determine an optimum working concentration for use in particular application.

Validity: 3 years.

**Hazard/Biohazard:** This material contains Proclin 300 as a preservation of the material, including its distribution. The end user assumes all responsibility for care, custody and control of the material, including its distribution. There is no known biohazard associated with this proclination.

**Remark:** This product is intended for research and manufacturing uses only. It is a diagnostic device. The user assumes all responsibility for care, custod and control of the material, including its disposal, in accordance with all regulations.

