

## Product Data Report

**Product Name:** Human IgE Antibody

**Catalog:** EE000502

**Source:** HEK 293 Cell

**Lot:** 211125

**Mfg Date:** 2021-11-25

**Exp Date:** 2024-11-24

**Con:** 5.2mg/ml

**Purity:** >95%, by SDS-PAGE, Affinity chromatography

**Buffer:** 1\*PBS

**Preservative:** 0.1% Proclin300

**Store:** Store at  $-20\pm 5$  °C. Repeated opening and closing of the bottle should be avoided. Antibodies are stable until the expiry date in an unopened vial.

**Application:** LF/IFA/ELISA/CMIA/WB

**Usage:** For application positive control. Each laboratory should determine an optimum working concentration for use in particular application.

**Validity:** 3years.

**Hazard/Biohazard:** This material contains Proclin 300 as a preservative. Although the amount is 0.1%, this substance is poisonous and hazardous and it should be handled by trained staff only. Take appropriate precaution. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations. There is no known biohazard associated with this product.

**Remark:** This product is intended for research and manufacturing uses only. It is a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

